

Шприц инъекционный B.Braun Injekt-F Luer Solo

ФСЗ 2012/12690

инъекционный, двухкомпонентный, объем 1 мл, с
иглой / без иглы, соединение Луер, положение канюли
центральное

B.Braun Melsungen AG, Германия



Характеристики

Тип	инъекционный
Исполнение	двухкомпонентный
Объем	1 мл
Игла	с иглой / без иглы
Соединение	Луер
Положение канюли	центральное

Дополнительные характеристики

Зеленый поршень удобен при определении объема наполнения.
Черная контрастная **градуировка 0.01 мл**.
Абсолютно прозрачный цилиндр шприца для контроля дозировки.
Исполнение **с иглой и без иглы**.

Injekt®-F Luer Solo, manufactured according to ISO 7886-1

Size	Cone	Graduation	Total Graduated Capacity	Units Per Box	Code No. (REF)
1 ml	Centric	0.01 ml	1.0 ml	100 pcs.	9166017V

Injekt®-F Luer Duo, manufactured according to ISO 7886-1 and ISO 7864

Size	Cone	Graduation	Total Graduated Capacity	Gauge x Length (Needle)	Ø x Length (Needle)	Units Per Box	Code No. (REF)
1 ml	Centric	0.01 ml	1.0 ml	25 G x 5/8 "	0.5 x 16 mm	100 pcs.	9166033V

В упаковке - 100 шт.